



# CLINICAL PATHOLOGY LABORATORIES

A Sonic Healthcare Clinical Laboratory

## Tube Guide Order of Draw

### Special Instructions

- Two patient identifiers must be on all specimens (full name, date of birth, ID#, etc.)
- The sequence below is the order the evacuated tubes must be filled to prevent specimen or additive cross contamination

For specimen requirements including storage and transport temperatures, consult our Test Directory on our website - [www.cpllabs.com](http://www.cpllabs.com)



Collection Containers	Bacterial Culture Media		Sodium Citrate	No Additives	Serum Separator Gel		Heparin	Potassium EDTA	Potassium EDTA Trace-Element Free	Sodium Fluoride or Potassium Oxalate	ACD Solution A or B	No Additive, Trace-Element Free
Additive	Bacterial Culture Media		Sodium Citrate	No Additives	Serum Separator Gel		Heparin	Potassium EDTA	Potassium EDTA Trace-Element Free	Sodium Fluoride or Potassium Oxalate	ACD Solution A or B	No Additive, Trace-Element Free
Notes	Invert 4 to 5 times		Invert 3 to 4 times tube must be filled completely	Invert 5 times blood should clot 60 minutes but no longer than 2 hours prior to separation	Invert 5 times blood should clot for 30 minutes but no longer than 2 hours prior to separation		Invert 8 to 10 times	Invert 8 to 10 times	Invert 8 to 10 times	Invert 8 to 10 times do not confuse with urine culture tube	Invert 8 to 10 times do not confuse with urine culture tube	Invert 8 to 10 times
Common Tests	<b>Order Code Name</b> 6013 Culture, Blood	<b>Order Code Name</b> 1425 Prothrombin Time (PT) 1430 PTT	<b>Order Code Name</b> 3061 Phenobarbital 3082 Theophylline 3025 Carbamazepine (Tegretol)	<b>Order Code Name</b> 142 Basic Metabolic Panel 9179 Comprehensive Metabolic Panel 115 Electrolyte Panel 9325 Hepatitis Panel, Acute 9175 Hepatic Function Panel 173 Lipid Panel 9324 Renal Function Panel  <b>Order Code Name</b> 3550 ANA (Anti-Nuclear Ab w/ Reflex Titer) 2075 CK, Total 3545 C-Reactive Protein (CRP) 2214 Creatinine 2675 Estradiol 2090 Ferritin 2700 Follicle Stim Hormone (FSH) 2217 Glucose 2713 HCG, Quantitative 2737 Hepatitis B Surface Ab	<b>Order Code Name</b> 2739 Hepatitis B Surface Ag 4675 Hepatitis C Ab 3540 HIV 1/2 4th Gen, Rfx Conf 4288 Homocysteine 4592 Herpes Simplex 1/2 IgG 2760 Insulin 2776 Luteinizing Hormone 3510 Mononucleosis Screen 2790 Progesterone 2800 Prolactin 2606 PSA, Total 3505 Rheumatoid Factor, Quant 3500 RPR 4600 Rubella Ab Screen 4273 T3, Free 2823 T4, Free 2819 T4, Total 2830 Testosterone 4513 Thyroid Peroxidase Ab 2835 TSH 2233 Uric Acid 2840 Vitamin B-12 4958 Vitamin D, 25 OH	<b>Lithium Heparin</b> <i>For Lithium Heparin tests, consult our Test Directory on our website: <a href="http://www.cpllabs.com">www.cpllabs.com</a></i>  <b>Sodium Heparin</b> <b>Order Code Name</b> 5207 Chromosome Analysis, Bone Marrow 5240 Chromosome Analysis, Peripheral Blood 8180 Leukemia/ Lymphoma Phenotyping, Peripheral Blood 8180 Leukemia/ Lymphoma Phenotyping, Bone Marrow 5320 Lymphocyte Antigen/Mitogen 4956 Vitamin B-6	<b>Order Code Name</b> 3810 Antibody Screen 3800 Blood Group (ABO) and RH Type 1000 CBC w/ Auto Diff w/ Platelets 1041 CBC w/out Diff w/ Platelets 2708 Hemoglobin A1c 4864 HLA-B27 8180 Leukemia/ Lymphoma Phenotyping, Peripheral Blood 1055 Sedimentation Rate (ESR) 4952 Vitamin B1	<b>Order Code Name</b> 4203 Arsenic 5641 Chromium, RBC 5696 Copper, RBC 4239 Lead 5698 Manganese, RBC 4245 Mercury 5699 Selenium, RBC 5632 Zinc, RBC	<b>Order Code Name</b> 3040 Ethanol, Blood 2217 Glucose 2119 Lactic Acid, Plasma	<b>Solution A</b> <b>Order Code Name</b> 4999 Hexosaminidase A/Total  <b>Solution B</b> <b>Order Code Name</b> 3750 HLA Class I Typing Disease Assoc 3876 HLA Class II Typing Disease Assoc 5005 Narcolepsy Evaluation  <b>Solution A or B</b> <b>Order Code Name</b> 5855 Alpha-Galactosidase 4972 Arylsulfatase A 5706 Beta-Glucosidase	<b>Order Code Name</b> 4265 Chromium 4219 Copper 5705 Iodine 4243 Manganese 4247 Molybdenum 4257 Selenium 4261 Zinc	

Clinical Pathology Laboratories | 9200 Wall Street, Austin, TX 78754 | [www.cpllabs.com](http://www.cpllabs.com) | P: 512.339.1275